

Shine-mini Portable Coagulation Analyzer

Create excellent quality with ingenuity



Shine-mini Portable Coagulation Analyzer

Create excellent quality with ingenuity

PRODUCT FEATURES

A DROP OF BLOOD

Blood coagulation function can be detected via ultramicro-detection of only 20ul blood samples, which meet the needs of infants and patients with difficulties in blood collection.


D-LAG

ONE-MINUTE TEST

Take PT for example, with an efficiency of 60 tests per hour, every detection can be done in 1 minute after loading the sample, thus to shorten the waiting time.

ONE-KEY TEST

The intelligent software system is equipped with a 3.5-inch color touch panel. The detection can be initiated by pressing one button.



Item	Sample volume	Measuring Range	Clinical Reference	200 ACT 150 y = 1 0118y - 3 9005 R ² = 0 9903
PT	20ul	7-90.5	10-14S	50
INR	20ul	0.5-9	0.7-1.3	e 50 100 150 200
ARTT	20ul	18-120.5	22-38.5	36 PT
Π	20ul	10-905	14-20S	20 y = 0.9876w - 0.1792 R ³ = 0.9971
FIB	20ul	1.5-6g/L	2-4g/L	10
ACT	20ul	50-500S	80-1405	0 9 3 10 15 20 25



MEDICAL AND LAB SUPPLIER

Komp. Ruko Sunter Paradise Tahap II Blok A4 No.43 Sunter Agung, Tanjung Priok - Jakarta Utara 14350 Telp.: (021) 6530 8478 - 6583 2509, HP: 0813 8892 7118 Fax.: (021) 6530 8479 E-mail: indolab@cbn.net.id